

T·SENSOR

T-Sensor is the latest product for digital radiography, designed with round corner and smooth focusing on high resolution, durability and patient comfort.



USB Type



Technical Specification



Model Type	USB 1.5
Scintillator type	CsI(+FOP)
Active Area	20(H)*30(V)mm
Active Pixels	1000*1500px
Pixel Pitch	20 _μ m
Frame rate	0.9 frames/s
Resolution	20Lp/mm
A/D Conversion	14bits
Interface	Direct USB
Energy Range	50~100kVp
Dimension	27*39mm



Make the perfect connection



Image in high resolution



Ensure product quality



Ergonomic design