

A person wearing a white lab coat is holding a white, handheld medical device. The device has a cable attached to one end and a small blue circular button or indicator on the side. The person's hand is wearing a white glove. The background is a blurred clinical setting.

PANDA smart

Smaller and Lighter

221x27x25mm

Only 138g



Parameters

Faster and More Accurate

1 min Ultra-fast Scan

7 μm Ultra-high Accuracy

PRODUCT INFORMATION

Technology	Staring Scanning with Coded Structure Light
Anti Fog System	Intelligent Heating
Accuracy	≤7μm
Precision	≤5μm
Handpiece Weight	138g (192g with cable)
Handpiece Dimension	221x27x25mm
Handpiece Enclosure	Aviation Aluminum Alloy
Scan Field	18x16mm
Scan Depth	0-20mm
True Color	Full HD
Tips	x6 (Autoclavable)
Tip Autoclave Cycle	50-100 times
Scanning Control	Button
Motion Sensing	Support
Cooling	Turbo Fan
Calibrator	Automatic
Data Interface	USB3.0
Cable Length (m)	1.85m

AFTERSALES

Standard Warranty	1 year
Extended Warranty	2-3 years available

SOFTWARE

Types of Scanning	Orthodontics, Implant, Restoration
Image Store	Yes
3D Data Preview	Yes
Output Files	STL, PLY, PTY (encrypted for user)
Cloud Service	Yes
Cloud Service Price	5GB/free
Direct Compatibility with Dental CAD	Exocad*
Main Image Processing Tools	Trim/Edit/Recover/Lock, etc.
Margin Line Drawing	Yes
Occlusal Match	Yes
Undercut	Yes
Data Comparison	Yes
AI Technology	Auto Tissue Remove, Overlap Correct, etc.
Software License	Free
Software Upgrade	Free
Transfer Data with API	Available

*All names, logos, brands and other trademarks featured or referred to within this brochure are the property of their respective trademark holders.

PC REQUIREMENTS

CPU	Intel i7-8700(3.2GHz) and above
Memory	DDR2400 16G and above
Hard Disk	Solid state drive 120G and above
GPU	NVIDIA GTX1060 (6G) and above
Computer Operating System	Windows10*64 bit
Display	Resolution 1920×1080 and above



Plug and Play

No need for multiple cables and power adapter sockets! Connect to PC with just one cable!





Rainbow Colors





Extraordinary Experience

 Yellow

 Silver

 Purple

 Green

 Red

 White